

Analytical Standard Record

P6F1048

Description:	AR16602##	Expires:	06/03/2014
Standard Type:	Calibration Standard	Prepared:	05/08/2011
Solvent:	NA	Prepared By:	Restek (Lot#: A033265)
Final Volume (mls):	1	Department:	ORGANICS
Vials:	1	Last Edit:	05/08/2011 00:00 by DBA

Analyte	Parent	CAS Number	Concentration	Units
Aroclor-1016	P6F1042	12674-11-2	0.2	ug/mL
Aroclor-1016 (1)	P6F1042	NA	0.2	ug/mL
Aroclor-1016 (1) [2C]	P6F1042	NA	0.2	ug/mL
Aroclor-1016 (2)	P6F1042	NA	0.2	ug/mL
Aroclor-1016 (2) [2C]	P6F1042	NA	0.2	ug/mL
Aroclor-1016 (3)	P6F1042	NA	0.2	ug/mL
Aroclor-1016 (3) [2C]	P6F1042	NA	0.2	ug/mL
Aroclor-1016 [2C]	P6F1042	12674-11-2	0.2	ug/mL
Aroclor-1260	P6F1042	11096-82-5	0.2	ug/mL
Aroclor-1260 (1)	P6F1042	NA	0.2	ug/mL
Aroclor-1260 (1) [2C]	P6F1042	NA	0.2	ug/mL
Aroclor-1260 (2)	P6F1042	NA	0.2	ug/mL
Aroclor-1260 (2) [2C]	P6F1042	NA	0.2	ug/mL
Aroclor-1260 (3)	P6F1042	NA	0.2	ug/mL
Aroclor-1260 (3) [2C]	P6F1042	NA	0.2	ug/mL
Aroclor-1260 [2C]	P6F1042	11096-82-5	0.2	ug/mL
Decachlorobiphenyl	P6F1042	2051-24-3	0.02	ug/mL
Decachlorobiphenyl [2C]	P6F1042	2051-24-3	0.02	ug/mL
Tetrachloro-m-xylene	P6F1042	877-09-8	0.01	ug/mL
Tetrachloro-m-xylene [2C]	P6F1042	877-09-8	0.01	ug/mL

Parent Standards used:

Standard	Description	Prepared	Prepared By	Lot Nbr	Expires	Last Edit	(mls)
P6F1042	AR1660 Stock	05/08/2011	Restek	A033625	06/03/2014	05/08/2011 00:00 by DBA	0.02