

PREPARATION BENCH SHEET

Print Date/Time: 11/03/2014 12:48 pm

Metals

B4I0004

Matrix: Water

Prepared using: Metals - EPA 200.7

Lab Number	Prepared	Initial (mL)	Final (mL)	Spike ID	Source ID	ul Spike	Comments
B4I0004-BLK1	10/24/2014 15:05	50	50				
B4I0004-BS1	10/24/2014 15:05	50	50	PI41501		500	
B4I0004-MS1	10/24/2014 15:05	50	50	PI41501	P409002-01	500	
B4I0004-MSD1	10/24/2014 15:05	50	50	PI41501	P409002-01	500	
P409002-01	10/24/2014 15:05	50	50	Client: Central Copper and Wire			Added for BatchQC in: B4I0004 BatchQC Analytes: Al, As, Ba, Ca, Cd, Co, Cr, Cu, Pb, Se, Zn
P409002-02	10/24/2014 15:05	50	50	Client: Central Copper and Wire			Analytes: As, Cd, Cr, Cu, Pb, Se, Zn
P409002-03	10/24/2014 15:05	50	50	Client: Central Copper and Wire			Analytes: As, Cd, Cr, Cu, Pb, Se, Zn
P409003-01	10/24/2014 15:05	50	50	Client: Western Remediation			Analytes: Al, Ba, Ca, Cd, Co, Cr
P409003-02	10/24/2014 15:05	50	50	Client: Western Remediation			Analytes: Al, Ba, Ca, Cd, Co, Cr

Batch Comments:

This batch has element names for the analytes ...

Spiking Witnessed By _____

Date _____

Preparation Reviewed By _____

Date _____

Extracts Received By _____

Date _____