

Data Review Report

(Continued)

Printed: 11/06/2014 1:56 pm

Analysis	Analyte	Initial Result	Result	Units	True Value	Rcvry	RPD
Lab ID: 9011801-BLK1 Blank							
Ag Total ICP 200.7	Silver	0.44	0.443	ug/L			
Cd Total ICP 200.7	Cadmium	0.00	0.00	ug/L			
Cr Total ICP 200.7	Chromium	0.00	0.00	ug/L			
Lab ID: 9011801-BS1 LCS							
Ag Total ICP 200.7	Silver	19.60	19.6	ug/L	20.0	98.0	
Cd Total ICP 200.7	Cadmium	21.20	21.2	ug/L	20.0	106	
Cr Total ICP 200.7	Chromium	20.00	20.0	ug/L	20.0	100	
Lab ID: 9011801-MS1 Matrix Spike							
Ag Total ICP 200.7	Silver	19.60	19.6	ug/L	20.0	98.0	
Cd Total ICP 200.7	Cadmium	59.45	59.4	ug/L	20.0	99.0	
Cr Total ICP 200.7	Chromium	44.60	44.6	ug/L	20.0	96.0	
Lab ID: 9011801-MSD1 Matrix Spike Dup							
Ag Total ICP 200.7	Silver	19.00	19.0	ug/L	20.0	95.0	3.11
Cd Total ICP 200.7	Cadmium	60.85	60.8	ug/L	20.0	106	2.33
Cr Total ICP 200.7	Chromium	45.20	45.2	ug/L	20.0	99.0	1.34