



**1664 HEM FOG in Water (EPA 1664)**

**Preservation:** Add HCl to pH<2; Store cool at 4°C

**Container:** Amber Glass - 1000 mL (HCl)

**Amount Required:** 1000

**Hold Time:** 10 days

| Analyte           | MDL | Reporting Limit | Surrogate %Rec | Duplicate RPD | ----Matrix Spike----<br>%Rec | RPD  | --Blank Spike / LCS--<br>%Rec | RPD  |
|-------------------|-----|-----------------|----------------|---------------|------------------------------|------|-------------------------------|------|
| Oil &Grease (HEM) | 5   | 10 mg/L         |                | 0.00          | 15.00-115.00                 | 0.00 | 15.00-15.00                   | 0.00 |